

Choose your obturation solutions



WaveOne® Gold
paper points



WaveOne® Gold
gutta-percha points



ThermoFil for
WaveOne® Gold

MAILLEFER INSTRUMENTS HOLDING S.A.R.L.
CHEMIN DU VERGER 3
CH-1338 BALLAIGUES
SWITZERLAND
www.dentsplymailefer.com

© 2017 Dentsply Sirona, Inc. Rx Only
ST8/ B EN WIGO TIP 000 / 01/2017 - updated 04/2017



WaveOne® Gold

Surf the canal
with confidence

+
WE
KNOW
ENDO.

 Dentsply
Sirona
Endodontics

One sequence for the majority of your cases

Scouting

ReadySteel®



010 02 - 25mm

Glide Path

WaveOne® Gold Glider



015 02 - 25mm
Variable Taper

Shaping

WaveOne® Gold



Size Primary,
025 07 - 25mm
Variable Taper

Product Line



015

2%



020

7%



025

7%



035

6%



045

5%